



# Recombinant Human Zinc finger protein 816 (ZNF816)

<b>Product Code</b>	CSB-EP610840HU-B
<b>Abbreviation</b>	ZNF816
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0VGE8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLREEATKKS KEKEPGMALP QGRLTFRDVA IEFSLLEWKC LNPAQRALYR AVMLENYRNL EFVDSSLKSM MEFSSTRHSI TGEVIHTGTL QRHKSHHIGD FCFPEMKKDI HHFEFQWQEV ERNGHEAPMT KIKKLTGSTD RSDHRHAGNK PIKDQLGLSF HSHLPELHMF QTKGKISNQL DKSIGASSAS ESQRISCRLK THISNKYGKN FLHSSFTQIQ EICMREKPCQ SNECGKAFNY SLLRRHHIT HSREREYKCD VCGKIFNQKQ YIVYHHRCHT GEKTYKNEC GKTFTQMSSL VCHRRLHTGE KPYKNECGK TFSEKSSLRC HRRLLHTGEKP YKCNECGKTF GRNSALVIHK AIHTGEKPYK CNECGKTF SQ KSSLQCHHIL HTGEKPYKCE ECDNVYIRRS HLERHRKIHT GEGSYKCKVC DKVFRSDSYL AEHQRVHTGE KPYKCNKCGR SFSRKSSLQY HHTLHTGEKP YTCNECGKVF SRRENLARHH RLHAGEKPYK CEECDKVFSR RSHLERHRR IHTGEKPYKCK VCDKAFRSDS CLANHTRVHT GEKPYKCNKC AKVFNQKGIL AQHQRVHTGE KPYKNECGK VFNQKASLAK HQRVHTAEKP YKCNECGKAF TGQSTLIHHQ AIHGCRET LQ M
<b>Source</b>	E.coli
<b>Target Names</b>	ZNF816
<b>Protein Names</b>	Recommended name: Zinc finger protein 816
<b>Expression Region</b>	1-651
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.