



Recombinant *Oryza sativa* subsp. *japonica* Nicotianamine synthase 3 (NAS3)

Product Code	CSB-EP610915OFG
Abbreviation	NAS3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0D3F2
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MTVEVEAVTM AKKEEQPEEEE VIEKLVKIT GLAAAIGKLP SLSPSPEVNA LFTELVMTCI PPSSVDVEQL GAEAQDMRGR LIRLCADAEG HLEAHYSDVL AAHDNPLDHL ALFPYFNYYI QLAQLEYALL ARHLPAAPPP SRLAFLGSGP LPLSSLVLAA RHLPAASFHN YDICADANRR ASRLVRADRD LSARMAFHTS DVAHVTTDLA AYDVVFLAAL VGMAAEEKAR MVEHLGKHMA PGAALVVRSA HGARGFLYPV VDPEEIRRGF FDLAVHHPE GEVINSVIA RKPPVAAPAL EGGDAHAHGH GAVVSRPCQR CEMEARAHQK MEDMSAMEKL PSS
Source	<i>E.coli</i>
Target Names	NAS3
Protein Names	Recommended name: Nicotianamine synthase 3 EC= 2.5.1.43 Alternative name(s): S-adenosyl-L-methionine:S-adenosyl-L-methionine:S-adenosyl-methionine 3-amino-3-carboxypropyltransferase 3 Short name= OsNAS3
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.