



Recombinant Bovine Spermatogenesis-associated protein 2-like protein (SPATA2L)

Product Code	CSB-EP611182BO
Abbreviation	SPATA2L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0IIA6
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	≥85% (SDS-PAGE)
Sequence	MGSSSLSEDY RLCLERELRR GRAGVCGDPS LRAVLWHILV EDFDLHGALQ DDALALLTDG LWGRADLAPA LRGLARAFEL LELAAVHLYL LPWRKEFTTI KTFSGGYVHV LKGALSEDLL IQSFQKMGYV RRDARHLMVA ALPPARQLVQ VALGCFALRL ECEILGEVLA QLGTSVLP AE ELLQARRASV DVASCVAWLQ QRLAREEEPP PLPRRGSP TG CQARLDLYRD VQEDEGSDEA SLYGGPSPGP DSPTSELACQ PRFWEQSARL WGAGGGPWEP AEVSSPTSGA SEEEEPQPEA FSFLSLRREL LSRPGDLAPP HAPRSPEQAS PPIPEPPGY QMHTCLAPGA LPALCCDTCR QLHAAHCAAL PSCHPGHSLR TLRGNSQRRL WLQRAQVDAL LYDSPAAGP
Source	E.coli
Target Names	SPATA2L
Protein Names	Recommended name: Spermatogenesis-associated protein 2-like protein Short name= SPATA2-like protein
Expression Region	1-409
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.