



Recombinant *Xenopus tropicalis* Diablo homolog, mitochondrial (diablo)

Product Code	CSB-BP611189XBF
Abbreviation	diablo
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0IJ31
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus tropicalis</i> (Western clawed frog) (<i>Silurana tropicalis</i>)
Purity	>85% (SDS-PAGE)
Sequence	VL SVGVGTSVCA IPVGQRSEPT LSSESLIKRA VSLVADSSST FLSQTTYALV ESLTEYTTAV YTLISLQQKY TSLLDKMNSN EESAIWQVII GARVQMKQLK EQYLKYEESW QRAVSLSEMA AEAAYQSGAD QASVTVRNHI QIVQTQVQQA RNQAHIAEVQ LAASQTDEIK RTITEDKGNP PSGGSPRSSL SEEEEIPEAY LRED
Source	Baculovirus
Target Names	diablo
Protein Names	Recommended name: Diablo homolog, mitochondrial Alternative name(s): Direct IAP-binding protein with low pI Second mitochondria-derived activator of caspase Short name= Smac protein
Expression Region	39-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.