



Recombinant *Oryza sativa* subsp. *japonica* Ribulose biphosphate carboxylase small chain, chloroplastic (RBCS)

Product Code	CSB-YP611191OFG
Abbreviation	RBCS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0INY7
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MQV WPIEGIKKFE TLSYLPPLTV EDLLKQIEYL LRSKWVPCLE FSKVGFBVYRE NHRSPGYDYG RYWTMWKLPF FGCTDATQVL KELEEAKKAY PDAFVRIIGF DNVRQVQLIS FIAYKPPGCE ESGGN
Source	Yeast
Target Names	RBCS
Protein Names	Recommended name: Ribulose biphosphate carboxylase small chain, chloroplastic Short name= Ribulose biphosphate carboxylase small chain C Short name= RuBisCO small subunit Short name= RuBisCO small subunit C EC= 4.1.1.39
Expression Region	48-175
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.