



Recombinant *Oryza sativa* subsp. *japonica* Cyclin-D3-2 (CYCD3-2)

Product Code	CSB-EP611195OFG-B
Abbreviation	CYCD3-2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0J3H7
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAFATLFDL YCPEEHLDF HDTAADDLH LDLHLHQPPP PPPLDLDLP ALFHALRGKE DPLRPAADD GYGGVSAREA AVGWALRAVA RLGFSALTA LAVAYLDRCF LGGALRLGDR PWMARLAAVA CVALAAKVEE TRVPVLLDLQ LCAAERADPN EAYVFEDKTV RRMELLVLSA LGWRMHPVTP LSYLQPLLGT AHAARLHHCD TALLALMPDW RWPRHRPSAW AAAALLATAG WCGGGGGDDA ELLALIDAPK DEMAECAKII SEEA AAAAAG GIVIGGENKR KGAAGLYSAP ASPSGVIGAS ACFSCDSSSS SVDSLFAALE PPGRPIKRG AAATTADPLP ADEESRDWP PYAA
Source	E.coli
Target Names	CYCD3-2
Protein Names	Recommended name: Cyclin-D3-2 Alternative name(s): G1/S-specific cyclin-D3-2 Short name= CycD3;2
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.