



Recombinant *Oryza sativa* subsp. japonica Squamosa promoter-binding-like protein 2 (SPL2)

Product Code	CSB-EP611203OFG
Abbreviation	SPL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0JG11
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDWDAKMPSW DLGTVVGPSTG GGGGGGGGGG ALDLKLGAPT SWKTTTTVSA ASAAPAAVAP PPPPPASSSS SAAAAGKRAR AGQQQQAAMP ACSVEGCAAD LSKCVRDYHR RHKVCCEAHSK TAVVTVAGQQ QRFCQQCSRFR HLLGEFDEEK RSCRKRLDGH NKRRRKPQPD PLNPGNLFAN HHGAARFTSY PQIFSTAASM SPQETKWPAV VVKTEAADVF QEPYYHALHL NGAGAAAAAS IFHHGGNKAR KHHPFLTAD HGGGAAAASP LFGCQPFTIT PSSESRSSSS SRHSNGKMFA HDGGLDNCAL SLLSDNPTPT AQITIPQPLF AGGGQYGGGG GGDVSLTGLS YVRMAGKDTS ILAKSATTTA TTATPTTTTS AQLQYHGYYH HHVSADQGSS DAAIQALPFS SW
Source	E.coli
Target Names	SPL2
Protein Names	Recommended name: Squamosa promoter-binding-like protein 2
Expression Region	1-412
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.