



Recombinant Danio rerio Protein phosphatase 1 regulatory subunit 35 (ppp1r35)

Product Code	CSB-MP611274DIL
Abbreviation	ppp1r35
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0P427
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MEVCGEPLIR AAPEPLRDER IQTAELLCAD LDLSVSLTPE RPADRRRRQV RFNVDPLIT VNPEPRNNQP TARKPPNKDE HGVETDREQS RECDGQQTHE AELNTTLALR AELEEEAEQT FDAEKAVREK LQSSTLTKNH VNSKAAEGLN FPRSQQLYRA LVSVLSLRDQ LISQALQDRP ALAPPTASQN NKFSSPPPEG PDILQFYSPD KMLRETPLLP GDHIPLRPRR PVPRPAHTTF HLHRLHKLWE S
Source	Mammalian cell
Target Names	ppp1r35
Protein Names	Recommended name: Protein phosphatase 1 regulatory subunit 35
Expression Region	1-251
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.