



# Recombinant Bovine Mitochondrial-processing peptidase subunit alpha (PMPCA)

<b>Product Code</b>	CSB-BP611294BO
<b>Abbreviation</b>	PMPCA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0P5M8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SSGGAYP NIPLSSPLPG VPKPVFATVD GQEKFKTKVT TLDNGLRVAS QNKFGQFCTV GILINSGSRY EAKYLSGIAH FLEKLAFSST ERFDSKDEIL LTLEKHGGIC DCQTSRDTTM YAVSADSKGL DTVVGLLADV VLHPRLTDEE IEMARMAVQF ELEDLNMRPD PEPLLTEMVH EAAYRENTVG LHRFCAENV GKMDRDVLHA YLRNYYPDR MVLAVGVGEH AQLVECARCY LLGTCPAWGT GAAVHVDRSV AQYTGIVKL ERDMSNVSLG PTPFPELTHI MIGLESCSFL EGDFIPFAVL NMMMGGGGSF SAGGPGKGMF TRLYLNLVLR HHWMYNATSY HHSYEDTGLL CIHASADPRQ VREMVEIVTR EFVLMAGTVD VVELERAKTQ LTSMLMMNLE ARPVIFEDVG RQVLATRSRK LPHELICALIR DVKPEDIKRV ASKMLRGKPA VAALGDLSEL PAYEHVQAAL ASRDGRLPRV YRLFR
<b>Source</b>	Baculovirus
<b>Target Names</b>	PMPCA
<b>Protein Names</b>	Recommended name: Mitochondrial-processing peptidase subunit alpha EC=3.4.24.64 Alternative name(s): Alpha-MPP
<b>Expression Region</b>	34-525
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.