



Recombinant Bovine Mitochondrial-processing peptidase subunit alpha (PMPCA)

Product Code	CSB-MP611294BO
Abbreviation	PMPCA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0P5M8
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	SSGGAYP NIPLSSPLPG VPKPVFATVD GQEKFKTKVT TLDNGLRVAS QNKFGQFCTV GILINSGSRY EAKYLSGIAH FLEKLAFSST ERFDSKDEIL LTLEKHGGIC DCQTSRDTTM YAVSADSKGL DTVVGLLADV VLHPRLTDEE IEMARMAVQF ELEDLNMRPD PEPLLTEMVH EAAYRENTVG LHRFCAENV GKMDRDVLHA YLRNYYPDR MVLAVGVGEH AQLVECARCY LLGTCPAWGT GAAVHVDRSV AQYTGIVKL ERDMSNVSLG PTPFPELTHI MIGLESCSFL EGDFIPFAVL NMMMGGGGSF SAGGPGKGMF TRLYLNLVLR HHWMYNATSY HHSYEDTGLL CIHASADPRQ VREMVEIVTR EFVLMAGTVD VVELERAKTQ LTSMLMMNLE ARPVIFEDVG RQVLATRSRK LPHELICALIR DVKPEDIKRV ASKMLRGKPA VAALGDLSEL PAYEHVQAAL ASRDGRLPRV YRLFR
Source	Mammalian cell
Target Names	PMPCA
Protein Names	Recommended name: Mitochondrial-processing peptidase subunit alpha EC=3.4.24.64 Alternative name(s): Alpha-MPP
Expression Region	34-525
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.