



Recombinant Shigella flexneri serotype 5b Phosphatidylserine decarboxylase proenzyme (psd)

Product Code	CSB-BP611449SGG
Abbreviation	psd
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0SXB9
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri serotype 5b (strain 8401)
Purity	>85% (SDS-PAGE)
Sequence	MLNSFKLSLQ YILPKLWLTR LAGWGASKRA GWLTKLVIDL FVKYYKVDMK EAQKPDASY RTFNEFFVRP LRDEVRPIDT DPNVLVMPAD GVISQLGKIE EDKILQAKGH NYSLEALLAG NYLMADLFRN GTFVTTYLSP RDYHRVHMPC NGILREMIYV PGDLFSVNHL TAQNVNLFARNERVICLFD TFGPMAQIL VGATIVGSIE TVWAGTITPP REGIIKRWTW PAGENDGSVA LLKGQEMGRF KLG
Source	Baculovirus
Target Names	psd
Protein Names	Recommended name: Phosphatidylserine decarboxylase proenzyme EC= 4.1.1.65 Cleaved into the following 2 chains: 1. Phosphatidylserine decarboxylase alpha chain 2. Phosphatidylserine decarboxylase beta chain
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.