



Recombinant Escherichia coli O6:K15:H31 Siroheme synthase (cysG)

Product Code	CSB-BP611565EGY
Abbreviation	cysG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TC92
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MDHLPIFCQL RDRDCLIVGG GDVAERKARL LLDAGARLTV NALAFIPQFT AWADAGMLTL VEGPFDESLL DTCWLAI AAT DDDALNQRVS EAAESRRIFC NVVDAPKAAS FIMPSIIDRS PLMVAVSSGG TSPVLARLLR EKLESLLPLH LGQVAKYAGQ LRGRVKQQFA TMGERRRFWE KLFVNDRLAQ SLANNDQKAI TETTEQLINE PLDHRGEVVL VGAGPGDAGL LTLKGLQQIQ QADV VVYDRL VSDDIMNLVR RDADRVFVGK RAGYHCVPQE EINQILLREA QKGKRVVRLK GGDPFIFGRG GEELETLCNA GIPFSVVPGI TAASGCSAYS GIPLTHRDYA QSVRLITGHL KTGGE LDWEN LAEKQTLVF YMGLNQAATI QKCLIEYGMP GEMPVAIVEN GTAVTQRVID GTLTQLGELA QQMNSPSLII IGRVVGLRDK LNWFSNH
Source	Baculovirus
Target Names	cysG
Protein Names	Recommended name: Siroheme synthase Including the following 3 domains: Uroporphyrinogen-III C-methyltransferase Short name= Urogen III methylase EC= 2.1.1.107 Alternative name(s): SUMT Uroporphyrinogen III methylase S
Expression Region	1-457
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.