



# Recombinant Escherichia coli O6:K15:H31 Multifunctional CCA protein (cca)

<b>Product Code</b>	CSB-MP611577EGY
<b>Abbreviation</b>	cca
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0TD50
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKIYLVGGAV RDALLGLPVK DRDWVVVGST PQEMLDAGYQ QVGRDFPVFL HPQTHEEYAL ARTERKSGSG YTGFTCYAAP DVTLEDDLKR RDLTINALAQ DDNGEIIDPY NGLGDLQNR LRVHSPAFGE DPLRVLRVAR FAARYAHLGF RIADETLALM REMTHAGELE HLTPERVWKE TESALTTRNP QVFFQVLRDC GALRVLFPEI DALFGVPAPA RWHPEIDTGI HTLMTLSMAA MLSPQVDVRF ATLCHDLGKG LTPPELWPRH HGHGPAGVKL VEQLCQRLRV PSEIRDLARL VAEFHDLIHT FPMLNPKTIV KLFDSIDAWR KPQRVEQLAL TSEADVGRGT GFESADYPQG RWLREAWEVA QSVPTKAVVE AGFKGVGIRE ELTRRRIAIV ASWKEQRCPK PD
<b>Source</b>	Mammalian cell
<b>Target Names</b>	cca
<b>Protein Names</b>	Recommended name: Multifunctional CCA protein Including the following 4 domains: CCA-adding enzyme EC= 2.7.7.72 Alternative name(s): CCA tRNA nucleotidyltransferase tRNA CCA-pyrophosphorylase tRNA adenylyl-/cytidylyl-transfe
<b>Expression Region</b>	1-412
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.