



Recombinant Escherichia coli O6:K15:H31 23S rRNA (uracil (1939)-C (5))-methyltransferase RlmD (rlmD)

Product Code	CSB-EP611581EGY-B
Abbreviation	rlmD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TE74
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	AQFYSAKRR TTTRQIITVS VNDLDSFGQG VARHNGKTLF IPGLLPQENA EVTVTEDKKQ YARAKVVRRL SDSPERETPR CPHFGVCGGC QQQHASVDLQ QRSKSAALAR LMKHEVSEVI ADVPWGYRRR ARLSLNYLPK TQQLQMGFRK AGSSDIVDVK QCPILVPQLE ALLPKVRACL GSLQAMRHLG HVELVQATSG TLMILRHTAP LSSADREKLE CFSHSEGLDL YLAPDSEILE TVSGEMPWYD SNGLRLTFSP RDFIQVNAGV NQKMVARALE WLDVEPEDRV LDLFCGMGNF TLPLATQAAS VVGVEGVPAL VEKGQQNARL NGLQNVTFYH ENLEEDVTKQ PWAKNGFDKV LLDPARAGAA GVMQQIikle PIRIVYVSCN PATLARDSEA LLKAGYTIAR LAMLDMPHPT GHLESMVLFS RVK
Source	E.coli
Target Names	rlmD
Protein Names	Recommended name: 23S rRNA (uracil(1939)-C(5))-methyltransferase RlmD EC= 2.1.1.190 Alternative name(s): 23S rRNA(m5U1939)-methyltransferase
Expression Region	2-433
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.