



Recombinant Escherichia coli O6:K15:H31 Glutamate--tRNA ligase (gltX)

Product Code	CSB-EP611589EGY
Abbreviation	gltX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TF64
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MKIKTRFAPS PTGYLHVGGG RTALYSWLFA RNHGGEFVLR IEDTDLERST PEAIEAIMDG MNWLNLEWDE GPYYQTKRFD RYNAVIDQML EEGTAYKCYC SKERLEALRE EQMAKGEKPR YDGRCRHSHE HHADDEPCVV RFANPQEGSV VFDDQIRGPI EFSNQELDDL IIRRTDGSPT YNFCVVVDDW DMEITHVIRG EDHINNTPRQ INILKALKAP VPVYAHVSMI NGDDGKKLSK RHGAVSVMQY RDDGYLPEAL LNYLVRLGWS HGDQEIFTRE EMIKYFALNA VSKSASAFNT DKLLWLNHHY INALPPEYVA THLQWHIEQE NIDTRNGPQL ADLVKLLGER CKTLKEMAQS CRYFYEDFAE FDADAAKKHL RPVARQPLEV VRDKLTAITD WTAENVHHAI QATADELEVG MGKVGMPPLRV AVTGAGQSPA LDVTVHAIGK TRSIERINKA LAFIAERENQ Q
Source	E.coli
Target Names	gltX
Protein Names	Recommended name: Glutamate--tRNA ligase EC= 6.1.1.17 Alternative name(s): Glutamyl-tRNA synthetase Short name= GluRS
Expression Region	1-471
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.