



Recombinant Escherichia coli O6:K15:H31 D-tagatose-1,6-bisphosphate aldolase subunit GatY (gatY)

Product Code	CSB-EP611597EGY-B
Abbreviation	gatY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TFZ6
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MYVVSTKQML NNAQRGGYAV PAFNIHNLET MQVVTETAAS MHAPVIIAGT PGTFTHVGTE NLMALVSAMA KQYHHPLAIH LDHHTKFDDI AQKVRSGVRS VMIDASHLPF AQNISRVKEV VDFCHRFDVS VEAELGQLGG QEDDVQVNEA DAFYTNPVQA REFAEATGID SLAVAIGTAH GMYARAPALD FSRLNIRQW VNLPLVLHGA SGLSTKDIQQ TIKLGICKIN VATELKNAFS QALKNYLTEH PEATDPRDYL QSAKSAMRDV VSKVIADCGC EGRA
Source	E.coli
Target Names	gatY
Protein Names	Recommended name: D-tagatose-1,6-bisphosphate aldolase subunit GatY Short name= TBPA Short name= TagBP aldolase EC= 4.1.2.40 Alternative name(s): D-tagatose-bisphosphate aldolase class II Tagatose-bisphosphate aldolase
Expression Region	1-284
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.