



# Recombinant Escherichia coli O6:K15:H31 HTH-type transcriptional repressor PurR (purR)

<b>Product Code</b>	CSB-BP611606EGY
<b>Abbreviation</b>	purR
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0THG9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MATIKDVAKR ANVSTTTVSH VINKTRFVAE ETRNAVWAAI KELHYSPSAV ARSLKVNHTK SIGLLATSSE AAYFAEIIIEA VEKNCQKGY TLILGNWNN LEKQRAYLSM MAQKRVDGLL VMCSEYPEPL LAMLEEYRHI PMVMDWGEE KADFTDAVID NAFEGGYMAG RYLIERGHRE IGVIPGLER NTGAGRLAGF MKAMEEAMIK VPESWIVQGD FEPESGYRAM QQILSQSHRP TAVFCGGDIM AMGALCAADE MGLRVPQDVS LIGYDNRNA RYFTPALTTI HQPKDSLGET AFNMLLDRIV NKREEPQSIE VHPRLIERS VADGPFRDYR R
<b>Source</b>	Baculovirus
<b>Target Names</b>	purR
<b>Protein Names</b>	Recommended name: HTH-type transcriptional repressor PurR Alternative name(s): Pur regulon repressor Purine nucleotide synthesis repressor
<b>Expression Region</b>	1-341
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.