



Recombinant Escherichia coli O6:K15:H31 Peptidyl-tRNA hydrolase (pth)

Product Code	CSB-YP611611EGY
Abbreviation	pth
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TIG6
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MTIKLIVGLA NPGA EYAATR HNAGAWFVDL LAERLRAPLR EEAKFFGYTS RVTLLGGEDVR LLVPTTFMNL SGKAVAAMAS FFRINPDEIL VAHDELDLPP GVAKFKLGGG HGGHNLKDI ISKLGNNPNF HRLRIGIGHP GDKNKVVG LVGKPPVSEQK LIDEAIDEAA RCTEMWFTDG LTKATNRLHA FKAQ
Source	Yeast
Target Names	pth
Protein Names	Recommended name: Peptidyl-tRNA hydrolase Short name= PTH EC= 3.1.1.29
Expression Region	1-194
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.