



Recombinant Escherichia coli O6:K15:H31 HTH-type transcriptional repressor AllR (allR)

Product Code	CSB-MP611629EGY
Abbreviation	allR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TKD5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MTEVRRRRGRP GQAEPVAQKG AQALERGIAl QYLEKSGGS SSVSDISLNL DLPLSTTFRL LKVLQAADFV YQDSQLGWWH IGLGVFNVGA AYIHNRDVLV VAGPFMRRLM LLSGETVNVA IRNGNEAVLI GQLECKSMVR MCAPLGSRLP LHASGAGKAL LYPLAEEELM SIILQTGLQQ FTPTTLVDMP TLLKDLEQAR ELGYTVDKKE HVVGLNCIAS AIYDDVGSVV AAISISGPSS RLTEDRFVSQ GELVRDTARD ISTALGLKAH P
Source	Mammalian cell
Target Names	allR
Protein Names	Recommended name: HTH-type transcriptional repressor AllR Alternative name(s): Negative regulator of allantoin and glyoxylate utilization operons
Expression Region	1-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.