



Recombinant Escherichia coli O6:K15:H31 Shikimate kinase 2 (aroL)

Product Code	CSB-BP611632EGY
Abbreviation	aroL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TKQ3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MTQPLFLIGP RGCGKTTVGM ALADSLNRRF VDTDQWLQSQ LNMTVAEIVE REEWAGFRAR ETAALEAVTA PSTVIATGGG IILTEFNRFH MQNNGIVVYL CAPVSVLVNR LQAAPEEDLR PTLTGKPLSE EVQEVLEERD ALYREVAHII IDATNEPSQV ISEIRSALAQ TINC
Source	Baculovirus
Target Names	aroL
Protein Names	Recommended name: Shikimate kinase 2 Short name= SK 2 EC= 2.7.1.71
Expression Region	1-174
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.