



# Recombinant Escherichia coli O6:K15:H31 Deoxyguanosinetriphosphate triphosphohydrolase (dgt)

<b>Product Code</b>	CSB-EP611635EGY-B
<b>Abbreviation</b>	dgt
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0TLH1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAQIDFRKKI NWHRRYRSPQ GVKTEHEILR IFESDRGR II NSPAIRRLQQ KTQVFPLERN AAVRTRLTHS MEVQQVGRYI AKEILSRLKE LKLEAYGLD ELTGPFESIV EMSCLMHDIG NPPFGHFGEA AINDWFRQRL YPEDAESQPL TDDRCSVAAL RLRDGEPLN ALRRKIRQDL CHFEGNAQGI RLVHTLMRMN LTWAQVGGIL KYTRPAWWRG ETPETHHYLM KKPGYYLSEE AYIARLRKEL NLALYSRFPL TWIMEAADDI SYCVADLEDA VEKRIFTVEQ LYHHLHEAWG QHEKGSLSL VVENAWEKSR SNSLSRSTED QFFMYLRVNT LNKLVPYAAQ RFIDNLPAIF AGTFNHALLE DASECSDLLK LYKNVAVKHV FSHPDVEQLE LQGYRVISGL LEIYRPLLNL PLSDFTELVE KERVKRFPIE TRLFHKLSTR HRLAYVEAVS KLPSDSPEFP LWEYYYRCRL LQDYISGMTD LYAWDEYRRL MAVEQ
<b>Source</b>	E.coli
<b>Target Names</b>	dgt
<b>Protein Names</b>	Recommended name: Deoxyguanosinetriphosphate triphosphohydrolase Short name= dGTP triphosphohydrolase Short name= dGTPase EC= 3.1.5.1
<b>Expression Region</b>	1-505
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.