



Recombinant Escherichia coli O6:K15:H31 Glutamate-1-semialdehyde 2,1-aminomutase (hemL)

Product Code	CSB-MP611636EGY
Abbreviation	hemL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TLH7
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MSKSENLISA ARELIPGGVN SPVRAFTGVG GTPLFIEKAD GAYLYDVDGK AYIDYVGSWG PMVLGHNHPA IRNAVIEAAE RGLSFGAPTE MEVKMAQLVT ELVPTMDMVR MVNSGTEATM SAIRLARGFT GRDKIIFEG CYHGHADCLL VKAGSGALTG QPNSPGVPA DFAKHTLTCT YNDLASVRAA FEQYPQEIAC IIVEPVAGNM NCVPLPEFL PGLRALCDEF GALLIIDEVM TGFRVALAGA QDYVGVEPDL TCLGKIIGGG MPVGAFGGRR DVMDALAPTG PVYQAGTSLG NPIAMAAGFA CLNEVAQPGV HETLDELTR LAEGLREAAE EAGIPLVVNH VGGMFGIFFT DAESVTCYQD VMACDVERFK RFFHMLLDEG VYLAPSAFEA GFMSVAHSME DINNTIDAAR RVFAKL
Source	Mammalian cell
Target Names	hemL
Protein Names	Recommended name: Glutamate-1-semialdehyde 2,1-aminomutase Short name= GSA EC= 5.4.3.8 Alternative name(s): Glutamate-1-semialdehyde aminotransferase Short name= GSA-AT
Expression Region	1-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.