



# Recombinant *Phaeosphaeria nodorum* Flap endonuclease 1 (FEN1)

<b>Product Code</b>	CSB-EP611709EUG
<b>Abbreviation</b>	FEN1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0UZR3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Phaeosphaeria nodorum</i> (strain SN15 / ATCC MYA-4574 / FGSC 10173) (Glume blotch fungus) ( <i>Parastagonospora nodorum</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGIKHLYQLI QEHSPPDAIKT GEIKNQFGRK VAIDASMSIY SFLIAVRSDG QQLMSETGET TSHLMGLFYR TLRMVDNGIK PLYVFDGAPP KLKSGELAKR FQRKSEAHAA AEEAKETGTA EDVEKFSRRT VRVTREHNEE CRLLKLMGI PFIAPTEAE AQCATLARGG KVYAAASEDM DTLTFNTPIIL LRHLTFSEQR KEPILEIHLD KVLEGLEMER EQPQFIDLCI LLGCDYLDPI KGIGPSTALK LIRDHKTLEG VVAHIQSDQK KLTLPEDWPF ADARQLFLEP DVCPADAPEC DFKWEAPDIE GLVKFLVEEK HFNEDRVRNG AAKLQKNMKS AQQSRLEGFF KPVERTPEEK ASLKRKADEK LSEK KKKQKE EAKAKKQAKS KPRTAG
<b>Source</b>	E.coli
<b>Target Names</b>	FEN1
<b>Protein Names</b>	Recommended name: Flap endonuclease 1 Short name= FEN-1 EC= 3.1.-.- Alternative name(s): Flap structure-specific endonuclease 1
<b>Expression Region</b>	1-396
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.