



# Recombinant Uncharacterized protein T19C3.5 (T19C3.5)

<b>Product Code</b>	CSB-YP611898CXY
<b>Abbreviation</b>	T19C3.5
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q10011
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	FNN ETLKVRLNVP FFKFFSKNGH HVVDEEIPKI TIPEIRLPFS SSIGHGVVRT DDLKIEKFSS PKIDFKLSKT GINWWTTGGA IKLSGRWHAK FTEVITIRDK GWLNAYATGI QMNISAAAYQ LDGQPQVKIG DCTVQIQKLD VEIGGSVLSW LVNLFETPFS KLVKKVINAQ ACSAAKGILI EEANRFLHSL PSHVDVGANF YVDYFLTQNP HATSEFTEFD LAADIVYGKS LCHPINIGNW TALLALNILY SSDASSTPGM LTSWISVSIP NCLIHSAHQN QLVKLLVSKD IPVVEPYLRT SCGFLGLCIG KFFKLRNEY PNNHVDLFFH TYSTPYFVMS EKEGVMLNMS LAVDLFINPY AKTKQNILAR LVVDTFSTVE PFLNHTRIHG RLQNSTITAR VDFSNIGDIP KAFLSAFSSI LSMTAREAVR SVLGVGIPIP SYDNVTLADS STIEVFDQFL RANIDFEYK
<b>Source</b>	Yeast
<b>Target Names</b>	T19C3.5
<b>Protein Names</b>	Recommended name: Uncharacterized protein T19C3.5
<b>Expression Region</b>	18-489
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.