



Recombinant Uncharacterized protein F52C9.3 (F52C9.3)

Product Code	CSB-BP611914CXY
Abbreviation	F52C9.3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10123
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MSKIIGKIGR LPKTLYEHHK KTIFFSFLGY LGADWVYRWD RNQGIRKEYA KIAVKYGETT VSPETRPKRV FVLVNVEGNS RGCDFQFNKN ALPLFHLAGV QVDVVKADNQ AQLEALAGAV DTQEADILYV VGGDGTIGTV VTGIFRNREK AQLPVGFYPG GYDNLWLKRM LPSVFENSDD VRHACETAMA VIEDQKKSIV AFELTTEGST LAPEYGLGDV SAGWFRQIED TRKKFWYFSM AKRRWAYFWE MLKRGPAPIE CHVEYEETCA GCEKCRPKPI IEAPQWRWWH VLTGTPKYKN NDGQKDYTGI INEKCGEKHE LDTHGAEFLI ENEQMSDYSQ IRFRMGGTDI GRFGVISDGF KRCSEGIVGR STDEKFGYTD FLANSVAFKI SALPSYIHRL YISSNATPKD AELTDRQITI RGTQKKLDIF LPTSIRLEL
Source	Baculovirus
Target Names	F52C9.3
Protein Names	Recommended name: Uncharacterized protein F52C9.3
Expression Region	1-439
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.