



Recombinant *Schizosaccharomyces pombe* Probable citrate synthase, mitochondrial (cit1)

| | |
|--------------------------|--|
| Product Code | CSB-EP611949SXV-B |
| Abbreviation | cit1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q10306 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | <pre> ASTR RLASSLLSQA SLRSRQLNPL FTSSYSTRSS SLKDRLAELI PEKQAEIKKF RAEHGQDVIG EVTINQMYGG ARGVRSLIWE GSVLDPNEGI RFRGYTIPEC QKLLPSSPNG KQPLPESLFW LLVTGEIPTL SQVQALSADW AARSQLPKFV EELIDRCPPT LHPMAQFSLA VTALEHDSAF AKAYERGMNK HDYWKYEYED CMDLIAKTVP IAGRIYRNLY RDGVVAPIQM DKDHSYNFAN VLGFANNEEF VELMRLYLTI HADHEGGNVS AHTGHLVGSA LSSPFLSMAA SLNGLAGPLH GLANQEVNLF LITMKKEIGD DLSEETIKSY LWKLLNSGRV VPGYGHAVLR KTDPRYTAQR EFALEHLPKD PMFQLVSRLY EIVPGVLTEH GKTKNPYPNV DSHSGVLLQY YGLKEQSFYT VLFQVSRVTLG VASQLIWDRD LGLPIERPXS FSTEALKKMV ETK </pre> |
| Source | E.coli |
| Target Names | cit1 |
| Protein Names | Recommended name: Probable citrate synthase, mitochondrial EC= 2.3.3.1 |
| Expression Region | 17-483 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.