



# Recombinant Schizosaccharomyces pombe ADP-ribosylation factor GTPase-activating protein glo3 (glo3)

<b>Product Code</b>	CSB-BP611960SXV
<b>Abbreviation</b>	glo3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q10367
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MTATKEESQK LLTSLRSQRD NKVCFDCGAK NPTWSSTTFG IYLCLDCSAA HRNMGVHISF VRSTVLDSWT YAQLRVMRVG GNENARNYFK RHGGVSLNLS KDCRLKYSSK TAKQYLEKLG SLAVEDEANY PDILDMDFLS NTHEGSSAAD TTNEDDDFFS AWDKASVKKS DDNLDDKTDL ASTSSSVVVE SGEKDEPVVV TEEKTMVSP SRPDSTSTTK SKTSSISSAR ARPIRASSRP TASKLGASRP QKLGIKKANA DIDFDEFEKA VLSSSESAPTK KPAAVASKES TVDTLVDNGV EEVKESTSTT VQGKPVKPV LKSAASAKSTK SDDSNLNANF ARLGFGQFAA ASNARAKAAA KARELKKNEV NAPTYARDHF ASQKSISSDQ YFGRGSFDPE AAAEAQERLS SFRDATAISS KSYFGEEEDE NEEGESSHRP DSAYLRDIAE TATEDIEAIK VAIHQGAEKL SDFIQKVGAR YNF
<b>Source</b>	Baculovirus
<b>Target Names</b>	glo3
<b>Protein Names</b>	Recommended name: ADP-ribosylation factor GTPase-activating protein glo3 Short name= ARF GAP glo3
<b>Expression Region</b>	1-483
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.