



Recombinant *Schizosaccharomyces pombe* Chromatin structure-remodeling complex subunit rsc58 (rsc58)

Product Code	CSB-MP611963SXV
Abbreviation	rsc58
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10412
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MNAESCSLIL SYISSAPAGN LVLNCTATTN KQTLEANLKS AVYQRYPEFL EDWTNVCVHI IHSYFPKDAC YWNVDKLYTS VNRITTLEQS KVVGKEIKNG FTFDAENSTA SLENDLQPQF RSHALFMVGS AGPLFSSTAR TSRLDSRLPD GGIIAKPVAL LPTPSVANSQ EYTLDKLSPP STAKPPASVI EFNPSLAKLP TVKYLQSGPF SSIAPYKNSS SSVIPDSSFH SVACYRASSH YKEAPVEKSI DIDIQNNLS LLEEDSWTSV PIQGELVELN KLLQHLQLLQ NQRITSHNVL SDEERQISVQ VQNLILKLAK DYDMSPEDFL MDDFTLSLTQ YGAFYRGTLP LSAQPLELPS QQLLRQSNA ALRSNSLSMN GSLSPSSTNV PLQSYRRTTK SRR
Source	Mammalian cell
Target Names	rsc58
Protein Names	Recommended name: Chromatin structure-remodeling complex subunit rsc58 Alternative name(s): RSC complex subunit rsc58
Expression Region	1-403
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.