



Recombinant *Schizosaccharomyces pombe* FAM203 family protein C1020.12c (SPAC26F1.12c)

Product Code	CSB-EP611978SXV
Abbreviation	SPAC26F1.12c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10498
Product Type	Recombinant Protein
Immunogen Species	<i>Schizosaccharomyces pombe</i> (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSELIELVGF LHDQNPQVRM LAVQHLLPYT ARNHPQFSIW FHNDFEPVKD LKALLKDKPQ IASQAVTALV NVSQNEKVRK VLMDDEFLQL IFSIVTNPLH GLADLSCMLL CNLAKEEDFA RILDMQVPLR EFSLSKNIID QLMDLFVKGT DHGINEYANF DFLANVFADM TRFERGRKYF TTLQEYDHVI PASKLVVFTE HKSLLRRTGV AAIKNISFD IPFQKVLMD EGINVLPYLL LPLAGPEELS EEDMDGMFDE LQLLPDDKKR EPDHFIMKTL VETLVLLTAT REGREHMRRR KVYPIIRELH LNVDEEIRE VCDQLVQMLV RDEAPEELEH IQDNPPDEDD VIVEVD
Source	E.coli
Target Names	hgh1
Protein Names	Recommended name: FAM203 family protein C1020.12c
Expression Region	1-356
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.