



# Recombinant Schizosaccharomyces pombe FAM203 family protein C1020.12c (SPAC26F1.12c)

<b>Product Code</b>	CSB-MP611978SXV
<b>Abbreviation</b>	SPAC26F1.12c
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q10498
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSELIELVGF LHDQNPQVRM LAVQHLLPYT ARNHPQFSIW FHNDFEPVKD LKALLKDKPQ IASQAVTALV NVSQNEKVRK VLMDDEFLLQL IFSIVTNPLH GLADLSCMLL CNLAKEEDFA RILDMQVPLR EFSLSKNIID QLMDLFVKGT DHGINEYANF DFLANVFADM TRFERGRKYF TTLQEYDHVI PASKLVVFTE HKSLLRRTGV AAIKNISFD IPFQKVLMD EGINVLPYLL LPLAGPEELS EEDMDGMFDE LQLLPDDKKR EPDHFIMKTL VETLVLLTAT REGREHMRRR KVYPIIRELH LNVDEEIRE VCDQLVQMLV RDEAPEELEH IQDNPPDEDD VIVEVD
<b>Source</b>	Mammalian cell
<b>Target Names</b>	hgh1
<b>Protein Names</b>	Recommended name: FAM203 family protein C1020.12c
<b>Expression Region</b>	1-356
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.