



Recombinant Putative zinc finger protein B0310.2 (B0310.2)

Product Code	CSB-BP612017CXY
Abbreviation	B0310.2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10938
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MEKSQPLSEF FSKQLVDPVS DALRSHTCSE QEMTLLESAL VHARDSLLEE IQRSKISEDV TQTVDCLLSR VDELLEETRN FAPISPLSSE APDILPAMIL FNSLCSSKSL DLEPIEIKVE AMENPVVSPF IFSSDDSVKG MPSRKRAKRS LEDTVQMLSK ENSVSPPPPS LPFVLPQIPV PIPFSNAALM QRWLGHPLTN PAYLNAMFQH QPAPKPSEEQ FAALMKITMH AANFMKTVPQ LPVSNHFTES DAEDIKIDVE SDEGEIEVSP SPSTGDITEN ESSSSSTGPM ISPNCGDGDC ALEKPFICMH NNCGRKFANK FLLKKHMFH TGLRPHTCPH CHKKFNRKDN LLRHKKTHSP TSAHLGPILP KNPFHGIPQL PVLHNIALTQ GFSHFHALKM SLETGATAVR SLS
Source	Baculovirus
Target Names	B0310.2
Protein Names	Recommended name: Putative zinc finger protein B0310.2
Expression Region	1-413
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.