



Recombinant *Oryza sativa* subsp. japonica Ninja-family protein Os03g0419100 (Os03g0419100, LOC_Os03g30570)

Product Code	CSB-MP612031OFG
Abbreviation	Os03g0419100, LOC_Os03g30570
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10JI3
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MASRDFLGVF GGGGGERRAA NGSGSAVGGE SDEIELSLGL SLGGRFGTDM SPDAKRARLA RSSSIASVCS VSAADGDPSP AAPLPLLRTS SLPTETEEER WRRREMQNRR RLEARRKRLR RRISVGSSSV PNKPGREDGG DGAVNRLQLR RSIGSQGSSS ANPQDQGPDG SAICQSTEAR SPSTSDDTNQ NSALPPTAST GKPLNGTVTQ QSPLRTLGLS TMRTSSTGDI GKIMMDMPMV SSKVEGPNGR KIDGFLYKYR KGEDVSIMCV CHGKFHSPAE FVKHAGGGDV SNPLRHIVVN PPSVFL
Source	Mammalian cell
Target Names	Os03g0419100
Protein Names	Recommended name: Ninja-family protein Os03g0419100
Expression Region	1-307
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.