



Recombinant *Oryza sativa* subsp. japonica Thioredoxin-like 1-2, chloroplastic (Os03g0326500, LOC_Os03g21000)

Product Code	CSB-MP612036OFG
Abbreviation	Os03g0326500, LOC_Os03g21000
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10M18
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAAAAA TSALRIGRRS PFLGRRLAVG PRRSRPVPRN LVAPVQMNLA FAKATKWWEK GLQPNMREVE SAQDLVDSL T NAGDNLVIVD FFSPGCGGCR ALHPKICQIA EQNPDLVFLQ VNYEEHKSMC YSLHVHVL PF FRFYRGAQGR LCSFSC TNAT IKKFRDALAK HKPDRCSLGP TRGLEESELL ALAANKDLQF NYTKKPELVP SGDAAAQEL DRGSTKLSPP AKPLVKQGSE ERSLVSSGR
Source	Mammalian cell
Target Names	Os03g0326500
Protein Names	Recommended name: Thioredoxin-like 1-2, chloroplastic Alternative name(s): Lilium-type thioredoxin 1-2
Expression Region	35-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.