



Recombinant *Oryza sativa* subsp. *japonica* Putative zinc finger CCCH domain-containing protein 21 (Os03g0301500, LOC_Os03g18950)

Product Code	CSB-BP612038OFG
Abbreviation	Os03g0301500, LOC_Os03g18950
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10MN8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MKLGISGGAG EAMEGELSFV SPARSSCF SF EGGGSGSPTW VSTVEALLRS PTSSVSDGGG GGGGGYN SPA RASSPLQKQI PYCRDAGDFS SLTWASTLEK PLESPSSCIS DGRGGGFGSP TSAFPPEKLL ISPPTCVSDN RGVGNVGGFP SLPWASSLER LLTSPSSCVS DSRGVGNADG FPSLPWASSL EKPLTSPSSC VSDGRSGGYS SPLGASAERE REVREAEMLL RAI AERYDDC FLRLRDAAAE LSDLHRERLR LAAENLHLSL LLEELESEQR KQASAMAPPK LEEDEAAQGG APKSISIRSP GYLSQKPPQG QARPQRLRVR ASQAMEISHP NCLIFVMGNQ CSPKEAAAAG DEEDEEDKGG GEVEVEAYRQ GAAKTELCKN WERGACPYGA RCRFAHGLQE LRPVIRHPRY KTLPCQMFAA ASGCPYGHRC HFRHSPLRAA AAESFCY
Source	Baculovirus
Target Names	Os03g0301500
Protein Names	Recommended name: Putative zinc finger CCCH domain-containing protein 21 Short name= OsC3H21
Expression Region	1-457
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.