



# Recombinant Schizosaccharomyces pombe 40 kDa peptidyl-prolyl cis-trans isomerase (wis2)

<b>Product Code</b>	CSB-BP612076SXV
<b>Abbreviation</b>	wis2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q11004
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSTYAYFKIS IDGKIQPTIY FELFDNVVPK TVKNFASLCN GFEKDGRCLT YKGSRFHRVI KNFMLQGGDF TRGNGTGGES IYGEKFEDEN FELKHDKPFL LSMANAGPNT NGSQFFITTV PTPHLDGKHV VFGKVIQGKS TVRTIENLET KNDDPVVPVV IEFCGTCTKD QIEAPKPDVT GDSLEEFDD YEGDKSETAI FKIASDLKGI ANKQFAQQNL DTAVAKWQKA LRYLMEYPVP NDDSKESPFD WKEYNALRYS IYANLALVAL KQNKPQEAIR NANIVIEASN STELEKQKAY YRLGCAQGLL KNFEESEKAL AKAGNDPAIS KKLAEIRQKK KDYKKRQKKA YAKMFQ
<b>Source</b>	Baculovirus
<b>Target Names</b>	wis2
<b>Protein Names</b>	Recommended name: 40 kDa peptidyl-prolyl cis-trans isomerase Short name= PPIase EC= 5.2.1.8 Alternative name(s): Cyclophilin-40 Short name= CYP-40 Rotamase
<b>Expression Region</b>	1-356
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.