



Recombinant Putative tyrosine-protein kinase C03B1.5 (C03B1.5)

Product Code	CSB-MP612104CXY
Abbreviation	C03B1.5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q11112
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MQSAQKSKTN GDTSLFLMKD IPQYWSPALK IGSGAFGEVM LVHSDRLDAR GEQILHNEFA LKVIRFPECD YHKIIREVNT HTLCRHHDNV LFIGSIYCQT FQQGLRVQMC LEYAALSDLS LLASSETFES HVAYICKNLI SALIHINLGL IVHGDLVKN ILMTHEGVVK LSDFGMANTF EQTRKSKNQS ILGTPGFIAP EIINLQGYDG KADLWGLGIL SLFLLKGENP FKKCTQFDLE SYKVSISESF YPDYSNYSTP LQEFMSSLKD YNPEFRSSAL HASSQTYLKT SCSQSEILNY YKRLRRTRGM DLPWPEEIHT FSVVPSLE
Source	Mammalian cell
Target Names	C03B1.5
Protein Names	Recommended name: Putative tyrosine-protein kinase C03B1.5 EC= 2.7.10.2
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.