



# Recombinant Protein let-756 (let-756)

<b>Product Code</b>	CSB-EP612114CXY-B
<b>Abbreviation</b>	let-756
<b>Storage</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
<b>Uniprot No.</b>	Q11184
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAVPAASSIV SYGGAATSNF LTTTPVPFLA GFYNSNFVTD RINSCAPYRV DRIRKQLQDE EENGYPADD RRRGALFCRS GTWLEMLPIE NPDDGSTRVK VHGTKEESSK FSIVEFVSV A MSLVSIRGVE TKNFICMDPS GKLYATPSSN YSTECVFLEE MMENYYNLYA SCAYGDRFNP WYIELRRSGK PRRGPNSKKR RKASHFLVVH HDLDRLRSPV PNGNDVTDLV VASLFHQPPS HPLFRQQTVT KPPNPHRISN LRAKVEMTNQ AEKQRLLEEK KRRREKKKRR REDRLRKEEQ IREARRQELK SLREEELRRR YQQQQQQQAS TQTRYNRPQN PANPYPTYRP LPTRSTVQSP RPAYNPYWQS PVTQAPHHNS HHHHHHHPRV SSSSDPQQRH QSQQHYLAQT VSNPNRQNVN YQRYP
<b>Source</b>	E.coli
<b>Target Names</b>	let-756
<b>Protein Names</b>	Recommended name: Protein let-756 Alternative name(s): Lethal protein 756
<b>Expression Region</b>	1-425
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>