



Recombinant *Saccharomyces cerevisiae* Protein SRL2 (SRL2)

Product Code	CSB-EP612269SVG-B
Abbreviation	SRL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12020
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSNFKNFTLN SFEDYYGKPS ETPKMEEEEKL EVTNVNASSS KKVHKSCKKST SKYDQKNVFR NSMTGIAQIL PTKPVKIIIEQ NIDFANPKSF DLLQSTHTIC FNKRINTTNT KLVNETHSS DIDNDILHVG APTDLGGNSN DEATRQLRK FRWSNNKEKS LCEKLTVIYW ALLLHTTKRA SKRRPILCHQ MIAEFFNRVY KEKSRVPITS RYIRDNLVAW VTQGKELHEK GWVGDAGTGD LQEQFNIATV KLYESAEDGR LSIGKDKPFR EENTGSDSLV RAEEDSTAVT NENGHISSEK NLKKDRRESI RNQILTLDLN DEDFFQNVMK VLSAIDPEL RQYVIVISEL VSMEMDDGKT VREKLRDVEL NINRLQVDIK EIKEMLVTLI NK
Source	E.coli
Target Names	SRL2
Protein Names	Recommended name: Protein SRL2 Alternative name(s): Suppressor of RAD53 null lethality protein 2
Expression Region	1-392
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.