



# Recombinant *Saccharomyces cerevisiae* Ribosome biogenesis protein YTM1 (YTM1)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP612270SVG   |
| <b>Abbreviation</b>      | YTM1  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q12024  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | <i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | MTEDKSQVKI RFFTREKDEL LHVQDTPMYA PISLKRYGLS EIVNHLLGSE<br>KPVPFDFLIE GELLRTSLHD YLTKKGLSSE ASLNVEYTRA ILPPSYLNSF<br>SNEDWVSSLD VGDGSKHIIS GSYDGIVRTW DLSGNVQKQY SGHSGPIRAV<br>KYISNTRLVS AGNDRTLRLW KTKNDDLKLT SQQQAQEDDD DEVNIEDGKT<br>LAILEGHKAP VVSIDVSDNS RILSASYDNS IGFWSTIYKE MTVVDPLEDI<br>NNPNNKISTA ARKRRKLTMK DGTIRRRAPL SLLESHTAPV EQVIFDSTDN<br>TVGYSVSQDH TIKTWDLVTA RCIDTRTTSY SLLSIAQLST LNLLACGSSA<br>RHITLH DPRV GASSKVTQQQ LIGHKNFVSS LDTCPENEYI LCSGSHDGTV<br>KVWDVRSTSP MYTITREDKS VQKGVNDKVF AVKWAKEVGI ISAGQDKKIQ<br>INKGDNIFKN |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | YTM1  |
| <b>Protein Names</b>     | Recommended name: Ribosome biogenesis protein YTM1 Alternative name(s):<br>Microtubule-associated protein YTM1  |
| <b>Expression Region</b> | 1-460   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life   |



of lyophilized form is 12 months at -20°C/-80°C.