



Recombinant *Saccharomyces cerevisiae* Transcription initiation factor TFIID subunit 10 (TAF10)

Product Code	CSB-MP612271SVG
Abbreviation	TAF10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12030
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MDFEEDYDAE FDDNQEGQLE TPFPSVAGAD DGDNDNDDSV AENMKKKQKR EAVVDDGSEN AFGIPEFTRK DKTLEEILEM MDSTPPIIPD AVIDYYLTKN GFNVADVVRVK RLLALATQKF VS DIAKDAYE YSRIRSSVAV SNANNSQARA RQLLQGQQQP GVQQISQQQH QQNEKTTASK VVLTVNDLSS AVAEYGLNIG RPDFYR
Source	Mammalian cell
Target Names	TAF10
Protein Names	Recommended name: Transcription initiation factor TFIID subunit 10 Alternative name(s): TAFII-23 Short name= TAFII23 TAFII-25 Short name= TAFII25 TBP-associated factor 10 TBP-associated factor 25 kDa p25
Expression Region	1-206
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.