



# Recombinant *Saccharomyces cerevisiae* Pre-mRNA-splicing factor CWC2 (CWC2)

<b>Product Code</b>	CSB-BP612274SVG
<b>Abbreviation</b>	CWC2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q12046
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTSWRDKSAK VQVKESELP S IPAQTGLTF NIWYNKWSQG FAGNTRFVSP FALQPQLHSG KTRGDNDGQL FFCLFFAKGM CCLGPKCEYL HHIPDEEDIG KLALRTEVLD CFGREKFADY REDMGGIGSF RKKNKTLTVG GIDGALNSKH LKPAQIESRI RFVFSRLGDI DRIRYVESKN CGFVKFKYQA NAEFAKEAMS NQTL LPSDK EWDDRREGTG LLVKWANEDP DPAAQKRLQE ELKLESLNMM VHLINNNTNS AGTEVNKNN ERLDRTFPEA SVDNVKKRLL PLDNGMESDD FIEK LKVKVKK NISRENISSK PSVGKLGGL LDYLSDED
<b>Source</b>	Baculovirus
<b>Target Names</b>	CWC2
<b>Protein Names</b>	Recommended name: Pre-mRNA-splicing factor CWC2 Alternative name(s): Complexed with CEF1 protein 2 PRP19-associated complex protein 40 Synthetic lethal with CLF1 protein 3
<b>Expression Region</b>	1-339
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.