



Recombinant *Saccharomyces cerevisiae* Pre-mRNA-splicing factor CWC2 (CWC2)

Product Code	CSB-YP612274SVG
Abbreviation	CWC2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12046
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MTSWRDKSAK VQVKESELP S IPAQTGLTF NIWYNKWSQG FAGNTRFVSP FALQPQLHSG KTRGDNDGQL FFCLFFAKGM CCLGPKCEYL HHIPDEEDIG KLALRTEVLD CFGREKFADY REDMGGIGSF RKKNKTYVVG GIDGALNSKH LKPAQIESRI RFVFSRLGDI DRIRYVESKN CGFVKFKYQA NAEFAKEAMS NQTL LPSDK EWDDRREGTG LLVKWANEDP DPAAQKRLQE ELKLESLNMM VHLINNNTNS AGTEVNNKNN ERLDRTFPEA SVDNVKKRLL PLDNGMESDD FIEKLLKVKK NISRENISSK PSVGKLGGPL LDYLSSDED
Source	Yeast
Target Names	CWC2
Protein Names	Recommended name: Pre-mRNA-splicing factor CWC2 Alternative name(s): Complexed with CEF1 protein 2 PRP19-associated complex protein 40 Synthetic lethal with CLF1 protein 3
Expression Region	1-339
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.