



Recombinant *Saccharomyces cerevisiae* Ribosome biogenesis protein NOP53 (NOP53)

Product Code	CSB-YP612279SVG
Abbreviation	NOP53
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12080
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MAPTNLTKKP SQYKQSSRKG KKAWRKNIDL SDVEQYMEKK IDHEITHGTS DITSLQNDAL FHVDVEGDEI LKNKLIKRRQ IKKVLKSKEI LDAVKTNSKI AALNHHKNSS GNPNKIQGVS KHELKLMAL AGRVHGESKI KNRVAKDGLV KTTAGDLWGE ESNSKKQKVK LPSGIKLDVE KKDQIPEELL KKSTTSWSTA SVRPSTLDIE PIAVKEFTEI PHAGKSYNPN NKAWSELINK EYKEEKARED ERIALEKYKE RIRHLMETLD DNEEEEESSN EEEEEEEEEEN ENENESTQCS GSDKEIKLSI NKPVKNNKKT KYQRNKAQRH EEKVKLQQL KELRQRVKDL EEVINSEETE ILSAIESDSN KVKKSKKNKK HCLGTYKYSVI DERLEIKFSD ELSDSLRKLK PEGNLLYDTV RKLQSSGKVE TRVPVRKGRK YKQKITEKWT HKDFK
Source	Yeast
Target Names	NOP53
Protein Names	Recommended name: Ribosome biogenesis protein NOP53 Alternative name(s): Nucleolar protein 53
Expression Region	1-455
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.