



Recombinant Schizosaccharomyces pombe Protein stn1 (stn1)

Product Code	CSB-BP612318SXV
Abbreviation	stn1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0E7J7
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MLENENQLTH QFPTLSRWNP MFISDVHKIS FHPHLQRYIG FWMGFPIRWI QIVGYIAAID IYEGKHVLTV DDCSGMVLRV VFIIQDDFSM SKRAISMSPG NVVCVFGKIN SFRSEVELIA QSFEELRDPN DEWKAWQKRM RYKKNLTKIS KNHHSIIRTP KKSYPKDHAK KELLKCIQRM CKLNVQAGFT IEELIYLKA KELHLPLVHI FNVENEIVEN CDGNHILALN FSLQTLQHG RIVRKSNSVY MLVTSKDIIR FVIPLMASGL LEAHKVQSIV RDSNPMFITL PLSAIAKHIC QFLLRTKGKW RQAKKYTWVR DNQFV
Source	Baculovirus
Target Names	stn1
Protein Names	Recommended name: Protein stn1
Expression Region	1-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.