



Recombinant *Oryza sativa* subsp. japonica Double-stranded RNA-binding protein 7 (DRB7)

Product Code	CSB-MP612568OFG
Abbreviation	DRB7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0IV63
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MDMPPTPLPP ETANTSPAPN GATAGIRVEN CYVFKSRLQE YAQKAGLQTP EYHTFKEGPS HEPVFKSTVV INNTSYDSL P GFFNRKAAEQ SAAEVALMEI VKSIPANANI PAVQETGLCK NLLQEYAQKM NYAIPSYICT KSASGLAPFI CTVEIGGIQY IGAARTKKD AEIKAARTAL LAIQGQSEGS ANGATKYIVV PGKRVGKEVE KMPIETPKPL KIKKGGFKKK WNKRFMKKD GQAVVEKDEA RVAGDAHSD VLMQPTVITQ EASCGTLFLQ PCEEAKRVEA EPPRDIEMVQ PDKENQHSDA ALVQPDDEAR VEQEPSRDIS VVPPNEEAIK VKQEPSIDAA ILQPKEEASS VKQEPFIDTA MLQACKEAGS VELGPARDTV ISQLNEQDRG VKQEPAGDTA VPQPDVDARV VKEESPRTEP NGEATNMKET PKNSAVCNSP ETKEFGDITA MGSDPPATNM SEE
Source	Mammalian cell
Target Names	DRB7
Protein Names	Recommended name: Double-stranded RNA-binding protein 7 Alternative name(s): dsRNA-binding protein 7
Expression Region	1-473
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.