



# Recombinant *Oryza sativa* subsp. japonica Double-stranded RNA-binding protein 7 (DRB7)

<b>Product Code</b>	CSB-YP612568OFG
<b>Abbreviation</b>	DRB7
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0IV63
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDMPPTPLPP ETANTSPAPN GATAGIRVEN CYVFKSRLQE YAQKAGLQTP EYHTFKEGPS HEPVFKSTVV INNTSYDSL P GFFNRKAAEQ SAAEVALMEI VKSIPANANI PAVQETGLCK NLLQEYAQKM NYAIPSYICT KSASGLAPFI CTVEIGGIQY IGAAARTKKD AEIKAARTAL LAIQGQSEGS ANGATKYIVV PGKRVGKEVE KMPIETPKPL KIKKGGFKKK WNKRFMKKD GQAVVEKDEA RVAGDAHSD VLMQPTVITQ EASCGLFLQ PCEEAKRVEA EPPRDIEMVQ PDKENQHSDA ALVQPDDEAR VEQEPSRDIS VVPPNEEAIK VKQEPSIDAA ILQPKEEASS VKQEPFIDTA MLQACKEAGS VELGPARDTV ISQLNEQDRG VKQEPAGDTA VPQPDVDARV VKEESPRTEP NGEATNMKET PKNSAVCNSP ETKEFGDITA MGSDPPATNM SEE
<b>Source</b>	Yeast
<b>Target Names</b>	DRB7
<b>Protein Names</b>	Recommended name: Double-stranded RNA-binding protein 7 Alternative name(s): dsRNA-binding protein 7
<b>Expression Region</b>	1-473
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.