



Recombinant *Oryza sativa* subsp. japonica CRS2-associated factor 2, mitochondrial (Os08g0188000, LOC_Os08g08860)

Product Code	CSB-EP612572OFG-B
Abbreviation	Os08g0188000, LOC_Os08g08860
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0J7J7
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	ATAAGE AIARRLNHHR APPFSDPDDD PPFTRLAERP PRAPSKKKKK EEEDQGGRIR PPEPASSDLP FDFRYSYSET DPAWRPIGFR EPTRFSPFGP GRLDRPWDGV AAAAARGEGA GAAATSREEV LGEPLAEEEV AQLVERYRHS DCSRQINLGK GGVTHNMIDD IHNHWKRAEA VRIKCLGVPT LDMDNICFHL EDKTGGKVIY RNINILILYR GRNYDPKQRP QIPLMLWKPL APIYPRLVQN VADGLTFEKT KELRNTGLNS SPLMKLTRNG VYVNVVDRVR EAFKTVEVVR LDCSHVGSSD CKKIGVKLRD LVPCVPLLFK DEQIILWRGK VKQENSVSLQ FSPEPS
Source	E.coli
Target Names	Os08g0188000
Protein Names	Recommended name: CRS2-associated factor 2, mitochondrial
Expression Region	15-366
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.