



# Recombinant *Oryza sativa* subsp. japonica Gibberellin 20 oxidase 2 (20ox2)

<b>Product Code</b>	CSB-YP612578OFG
<b>Abbreviation</b>	20ox2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0JH50
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVAEHPTPPQ PHQPPPMDST AGSGIAAPAA AAVCDLRMEP KIPEPFVWPN GDARPASAAE LDMPVVDVGV LRDGDAEGLR RAAAQVAAAC ATHGFFQVSE HGVDAALARA ALDGASDFFR LPLAEKRRAR RVPGTVSGYT SAHADRFASK LPWKETLSFG FHDRAAAPVV ADYFSSTLGP DFAPMGRVYQ KYCEEMKELS LTIMELLELS LGVERGYRE FFADSSSIMR CNYPPCPEP ERTLGTGPHC DPTALTILLQ DDVGGLEVLV DGEWRPVSPV PGAMVINIGD TFMALSNGRY KSCLHRAVVN QRRERRSLAF FLCPREDRVV RPPPSAATPQ HYPDFTWADL MRFTQRHYRA DTRTLDAFTR WLAPPAADAA ATAQVEAAS
<b>Source</b>	Yeast
<b>Target Names</b>	20ox2
<b>Protein Names</b>	Recommended name: Gibberellin 20 oxidase 2 EC= 1.14.11.- Alternative name(s): GA 20-oxidase 2 Gibberellin C-20 oxidase 2 Os20ox2 Protein semidwarf-1
<b>Expression Region</b>	1-389
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.