



# Recombinant Bovine Keratin, type I cuticular Ha5 (KRT35)

<b>Product Code</b>	CSB-MP612675BO
<b>Abbreviation</b>	KRT35
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0P5J7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MASKCLKASF SSGSLKVPGG AGGGSARVST IFSSSSCKLP SFSRGPRSFS ACSVGLGKSS CRAASCLPAL CLPSGGFATS YSMAGGWFGE GILTGNEKET MQFLNDRLAS YLEKCGSWG RTRSWRAAST SGVNSALPV PDYQSYFQTI EELQKKTLC T KSENARLVVQ IDNAKLAADD FR TKYETEVS MRQLVESDMN GLRRILDDLT LCKADLEA QV ESLKEELLCL KKNHEEEVNS LRCQLGDR LN VEVDAAPPVD LNRVLNEMRC QYETLVENNR REAEDWFNTQ TEELNQQVVS SSEQLQSYQA EI IELRRTVN ALEIELQAQH SMRDALESTL AETEAR YSSQ LAQMQLIGN VESQLAEIRC DLERQNQEYQ VLLDVRARLE CEINTYRGLL DSEDCKLPCN PCAPDHSPSK SCLPCLPAAS CGPGMARTTC SPRPICVPCP GSRF
<b>Source</b>	Mammalian cell
<b>Target Names</b>	KRT35
<b>Protein Names</b>	Recommended name: Keratin, type I cuticular Ha5 Alternative name(s): Keratin-35 Short name= K35
<b>Expression Region</b>	1-454
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This protein is a member of the keratin gene family. This type I hair keratin is an acidic protein which heterodimerizes with type II keratins to form hair and nails. The type I hair keratins are clustered in a region of chromosome 17q12-q21 and have the same direction of transcription.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.