



Recombinant Escherichia coli O6:K15:H31 HTH-type transcriptional repressor fabR (fabR)

Product Code	CSB-EP612939EGY-B
Abbreviation	fabR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TA93
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MGVRAQQKEK TRRSLVEAAF SQLSAERSFA SLSLREVARE AGIAPTSFYR HFRDVDELGL TMVDEGLML RQLMRQARQR IAKGGSVIRT SVSTFMEFIG NNPNAFRLLL RERSGTSAAF RAAVAREIQH FIAELADYLE LENHMPRAFT EAQAEAMVTI VFSAGAEALD VGVEQRRQLE ERLVLQLRMI SKGAYYWYRR EQEKTTIIPG NVKDE
Source	E.coli
Target Names	fabR
Protein Names	Recommended name: HTH-type transcriptional repressor fabR
Expression Region	1-215
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.